
# Hoefer SE300 miniVE

Mini-Vertikal-Gelelektrophorese-Einheit

# Inhalt

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## Important Information – English

- If this equipment is used in a manner not specified by Hoefer, Inc. the protection provided by the equipment may be impaired.
- This instrument is designed for indoor laboratory use only.
- Only accessories and parts approved or supplied by Hoefer, Inc. may be used for operating, maintaining, and servicing this product.
- Only use a power supply that is CE marked or safety certified by a nationally recognized testing laboratory.
- The safety lid must be in place before connecting the power supply leads to a power supply.
- Turn all power supply controls off and disconnect the power leads before removing the safety lid.
- Circulate only water or 50/50 water/ethylene glycol through the heat exchanger if so equipped. Do not connect the heat exchanger to a water tap or any coolant source where the water pressure is unregulated.
- Never introduce antifreeze or any organic solvent into any part of the instrument. Organic solvents will cause irreparable damage to the unit!
- Do not operate with buffer temperatures above the maximum specified technical specifications. Overheating will cause irreparable damage to the unit!

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## Bibliographie

### Polyacrylamid-Gelelektrophorese

- Adams, L.D. and Gallagher, S.R. 1992. Two-Dimensional Gel Electrophoresis Using the O'Farrell System. *Current Protocols in Molecular Biology*, 10.4.1–10.4.13.
- Gallagher, S.R., and Smith, J.A. 1991. Electrophoretic separation of proteins. *Current Protocols in Molecular Biology*. (F.A. Ausubel, R. Brent, R.E. Kingston, D.D. Moore, J.G. Seidman, J.A. Smith, and K. Struhl, eds.) 10.2.1–10.2.21.
- Laemmli, U.K. 1970. Cleavage of structural proteins during the assembly of the head of bacteriophage T4. *Nature*. 227, 680–685.
- Matsudaira, P.T. and Burgess, D.R. 1978. SDS microslab linear gradient polyacrylamide gel electrophoresis. *Anal. Biochem.* 87:386–396.
- Reisfeld, R.A., et al. 1962. Acidic buffer system for resolution of cationic proteins. *Nature*. 195:281.
- Sasse, J., and Gallagher, S.R. 1991. Staining proteins in gels. *Current Protocols in Molecular Biology*. (F.A. Ausubel, R. Brent, R.E. Kingston, D.D. Moore, J.G. Seidman, J.A. Smith, and K. Struhl, eds.) 10.6.1–10.6.8.
- Weber, K., and Osborn, M. 1969. The reliability of molecular weight determinators by dodecyl sulfate-polyacrylamide gel electrophoresis. *J. Biol. Chem.* 224:4406–4412.

### Blotting

- Gallagher, S., Winston, S.E., Fuller, S.A. and Hurrell, J.G.R. 1993. Immunoblotting and Immunodetection. in *Current Protocols in Molecular Biology*. 10.8.1–10.8.17. Greene Publishing and Wiley-Interscience, NY.
- Gershoni, J.M., and G.E. Palade. 1983. Protein Blotting: Principles and Applications. *Anal. Biochem.* 131:1–15.
- Harlow, E. and Lane, D. 1988. *Antibodies: A Laboratory Manual*, Cold Spring Harbor Press.
- Sambrook, J., et al. 1989. *Molecular Cloning: A Laboratory Manual*, Cold Spring Harbor Laboratory Press, B.23.
- Towbin, H., Staehelin, T., and Gordon, J. 1979. Electrophoretic transfer of proteins from polyacrylamide gels to nitrocellulose sheets: procedure and some applications. *Proc. Natl. Acad. Sci. USA*. 76:4350–4354.
-

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© 2012 Hoefer, Inc.
Warenzeichen von Hoefer, Inc. ist