
# MacroVue UV-25

Transilluminator

# Inhalt

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## Important Information – English

- If this equipment is used in a manner not specified by Hoefer, Inc. the protection provided by the equipment may be impaired.
- This instrument is designed for indoor laboratory use only.
- Only accessories and parts approved or supplied by Hoefer, Inc. may be used for operating, maintaining, and servicing this product.
- Warning! Because this instrument can develop sufficient voltage and current to produce a lethal shock, care must be exercised in its operation.
- This instrument is designed in accordance with the EN61010-1:2001 electrical safety standard. Nevertheless, it should be used only by properly trained operators. Read this entire manual before using the instrument and use only according to the instructions.
- The instrument must always be used with the earth lead of the power cord correctly grounded to earth at the mains outlet.
- Use only undamaged electrical wire and equipment specific for the voltages you will use. All equipment connected to high voltage should be in accordance with EN61010-1:2001.
- Keep the instrument as dry and clean as possible. Wipe regularly with a soft, damp cloth. Let the instrument dry completely before use.
- Do not operate the instrument in extreme humidity (above 80%). Avoid condensation by letting the unit equilibrate to ambient temperature when taking the instrument from a colder to a warmer environment.
- To permit sufficient cooling, ensure that the vents of the instrument are not covered.

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Khandjian, E. W. (1987) Optimized hybridization of DNA blotted and fixed to nitrocellulose and nylon membranes. *BioTechnology.* **5**, 165–167.

Khandjian, E. W. (1986) UV crosslinking of RNA to nylon membrane enhances hybridization signals. *Molec. Biol. Report.* **11**, 107–115.

Reed, K. C. and D. A. Mann. (1985) Rapid transfer of DNA from agarose gels to nylon membranes. *Nuc. Acid Res.* **13**, 7207–7221.

## DNA nicking

Vollrath, D. and R. W. Davis. (1987) Resolution of DNA molecules greater than 5 megabases by contour-clamped homogeneous electric fields. *Nuc. Acid Res.* **15**, 7865–7876.

## Rapid restriction site mapping

Whittaker, P. A. and E. M. Southern. (1986) Ultraviolet irradiation of DNA: A way of generating partial digests for rapid restriction site mapping. *Gene.* **41**, 129–134.

---

### Testing RecA function

Maniatis, T. (1989) Maintenance and Testing of Lysogens of Bacteriophage Lambda in *Molecular Cloning: A Laboratory Manual*, eds. Sambrook, Fritsch and Maniatis. 2.96.

### Ultraviolet sterilization

Reich, R. R., and Anderson, H. D. (1986) Sterilization of Membrane Filters with Ultraviolet Radiation. *Pharmaceutical Manufacturing*. Jan.

### Elimination of PCR contamination

Wilson, K. H. (1992) UV Absorption Complicates PCR Decontamination, *BioTechniques*. **13**, 208–210.

Sarkar, Gobinda, and Sommer, Steve S. (1991) Parameters Affecting Susceptibility of PCR Contamination to UV Inactivation, *BioTechniques*, **10**, 590–594.

Sarkar, Gobinda, and Sommer, Steve S. (1990) Shedding Light on PCR Contamination. *Nature*. **343**, 27.

Yamada, H. and Hieda, K. (1992) Wavelength Dependence (150–290nm) of the Formation of the Cyclobutane Dimer and the (6-4) Photoproduct of Thymine. *Photochemistry and Photobiology*. **55**, 541–548.

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© 2012 Hoefer, Inc.
Hoefer, Inc. bietet eine komplette technische Unterstützung für alle unsere Produkte.
Warenzeichen von Hoefer, Inc.